

WEARABLE DEVICE

Acesight



As all-in-one ease-to-use wearable smart glasses, Acesight VR is a fantastic option for staying connected – advanced visual technology ensures better reading and watching experiences.



Suitable for the following conditions:

- | | |
|------------------------|------------------------|
| ✓ Macular Degeneration | ✓ Nystagmus |
| ✓ Glaucoma | ✓ Esotropia |
| ✓ Diabetic Retinopathy | ✓ Retinitis Pigmentosa |
| ✓ Stargardt's Disease | ✓ And more |

01 Technical Specifications

Magnification	1X - 16X
Camera	48MP
Screen Resolution	Full HD 2340x1080
FOV	65°
Refresh Rate	60FPS
Weight of Headset	455g/1 pound
Weight of Controller	65g/0.14 pound
Dimensions of Headset	180x120x100mm/7x4.7x3.9in
Headset Charging Time	Less than 1.5 hours
Headset Battery Life	3.5-4 hours

Users' Reviews

Acesight VR is so crucial for my husband. He is a handicraftsman and he needs to stare at small details sometimes. Acesight VR helps him a lot so he can continue his work. We cannot thank you enough!

Dennis - March 2, 2022, USA

It's a fantastic gift for my Dad. He likes Acesight VR! It's perfect for him to wear it while watching TV shows. That was what he dreamed of. I am so glad that I selected the right thing for him.

Dorothy - November 13, 2022, Canada

