

# HANDHELD DEVICE

## Luna 6

### We Design It as You Need



Pocket-sized Luna 6 has a unique wireless charging system and built-in handle, offers sharp images even under fast movement. Users can easily carry it around and use it everywhere.

- Pocket-size and Lightweight
- Customized Wireless Charger
- Comfortable Built-in Handle



# 01 Technical Specifications

Magnification	2.5X - 19X
Screen	5.5-inch, Resolution 1280*720
Camera	13MP
Focus	Auto focus and refocus by short pressing the Mode button and Freeze button simultaneously
Color Modes	Full color and 10 high contrast color modes
Weight	Luna 6: 228g/0.5lbs; wireless charger: 92g/0.21lbs
Dimensions	160x80x11mm/6.3x3.1x0.4in
Battery Life	4 hours continuous use
Wired Charging Time	Less than 3 hours
Wireless Charging Time	About 4 hours

## 02 Users' Reviews

I bought Luna 6 from the clinic, and the service staff helped me a lot with the operation. I don't know how to express my gratitude to Zoomax. Thanks for making such a good product.

Tamara - April 5, 2022, Britain

I work in the library. I always need to check the book label. Luna 6 is the best device for me! It is light and slim so that I can bring it and use it anytime.

Kate - June 23, 2022, USA

