

# Snow Pad 2025

Our latest foldable and portable video magnifier. The split screens with near and distant viewing mode empowers low vision individuals to handle multi-tasks. Outstanding compatibility allows for seamless integration with your personal tablets.



**Ideal For Home, School, Work & Travel**  
**Compatible With Your Personal Tablet**

**ZOOMAX<sup>®</sup>**



Introducing the Zoomax Snow Pad, a next-generation portable digital magnifier that seamlessly integrates with your personal tablet. Designed for individuals with low vision, the Snow Pad transforms your tablet into a versatile visual aid, offering exceptional image quality for both distance and close-up viewing. Its features include split-screen functionality, OCR/text-to-speech capabilities, and a customized app tailored for ease of use. Whether you're at school, work, or enjoying leisure activities, the Snow Pad enhances your visual experience, promoting independence and accessibility in daily life.



Customized App



OCR/Text to Speech



File Edit



Share File



Split Screen



Magnification

- ① Download Snow Pad App to Your Tablet
- ② Connect Snow Pad with Your Tablet using The Supplied USB Cable
- ③ Power On The Snow Pad and Tap The Snow Pad App to Enjoy It

## Where to download the App?



Snow Pad iPadOS



Turn 270 degrees clockwise.  
Turn 60 degrees counterclockwise.



Up 15 degrees.  
Down 40 degrees.

|                          |  |
|--------------------------|--|
| Magnification            | 1X – 32X   |
| Display Modes            | Full color and 10 high contrast color modes                        |
| Build-in Distance Camera | Auto Focus, Optical Zoom 7X, 1920×1080, Full HD                    |
| Frame Rate               | 60FPS  |
| Folded Dimensions        | L*W*H/ 309*312*55mm  |
| Stand Height             | 278mm  |
| Weight                   | Approx. 2.4kgs (5.29lbs) without tablet                            |
| Battery                  | Approx. 5h   |
| Charging Time            | Approx. 4.5h   |
| Power Adapter            | Input 100–240V; Output 5V/2.5A                                     |
| OCR/Text-to-speech       | 20+ Languages  |
| Compatible W/Tablets     | iPads with iOS 16 or higher provide the best compatibility         |
| Viewing Modes            | 4 Options: Split Screen; Distance-viewing;<br>Near-viewing; Selfie |



[www.zoomax.com](http://www.zoomax.com)

[www.zoomaxusa.com](http://www.zoomaxusa.com)

[info@zoomax.com](mailto:info@zoomax.com)

